

# Quick reference guide

## SafeFast Classic vs SafeFast Premium



	<b>CLASSIC</b>	<b>PREMIUM</b>
<b>Class/Airflow</b>	Class II	Class II
<b>Product Type</b>	MBSC	MBSC
<b>Sizes</b>	3/4/5/6ft	3/4/5/6ft
<b>Double Fan</b>	Yes	Yes
<b>Single or Double HEPA Exhaust Filter</b>	Yes	Yes
<b>Recirculating or Ducted</b>	Yes	Yes
<b>Motorised Sliding and Hinged Front Window</b>	Yes Tempered glass	Yes Laminated glass
<b>Glass Side Windows</b>	Yes	Yes
<b>304 Stainless Steel Rear Wall</b>	No	Yes
<b>Powder Coated Rear Wall</b>	Yes	No
<b>304 Stainless Steel Work Top</b>	Yes	No
<b>316 Stainless Steel Work Top</b>	No	Yes
<b>LED Lights</b>	No Fluorescent tube	Yes
<b>DC Fans</b>	No AC	Yes
<b>HEPA Filters</b>	Normal H14	Low pressure drop H14
<b>1 x Electric Outlet (3/4ft models)</b>	Yes	No
<b>2 x Electric Outlet (5/6ft models)</b>	Yes	Yes
<b>1 x Gas Tap with solenoid Valve</b>	Optional	Optional
<b>1 x Vacuum Tap</b>	Optional	Optional